



New microLFA modules for on-line measurement of
nutrients in ship of convenience application

Dr. Ing. Luca Sanfilippo



User requirements for an on-line chemical analyzer in a Ferrybox

- Long term unattended autonomy
- Compactness
- Low reagents and sample consumption
- Low life-cycle cost
- Excellent reliability
- Easy interface with data-logger
- Low maintenance by non expert users





Ferrybox applications with Micromac-1000 analyzers



AWI – BAH
(Biologische Anstalt Helgoland)
2005



MUMM
(Belgica ship)
2011



Micromac-1000 features and limitations in Ferrybox systems

Features:


- Unattended long term use
- Low reagents consumption
- Sensitivity for sea water measurements
- Compactness and portability
- Modularity
- 12 Vdc power supply.

Limitations:

- Silicone based not sealed hydraulics
 - Limited internal space for reagents solutions
 - Hydraulics not directly visible to the user
 - Electronics not sealed.
-



The advance: MicroLFA module







MicroLFA module

- Very simple hydraulics directly accessible by operator
- Complete separation between hydraulic and sealed electronic
- Suitable to run all nutrients methods already developed and tested on field
- Lower reagents consumption: 60 μL for most of the reagents used
- Fast “plug-in” hydraulic connector allows easy deployment and reagents changeover on board
- Power supply: 12 Vdc, 3 W stand-by, 6 W analysis, max. 1 A
- Very low maintenance
- Suitable for long term deployment.

The micro Loop Flow Reactor



P- PO4 Analyzer – Seawater rev 0



Nutrients analytical methods using microLFA technology

- **AMMONIA (12 min.):**
OPA fluorimetric method, 6 ppb
- **NITRITE (5 min.):** NED-SAA, 1 ppb
- **NITRATE + NITRITE (12 min.):**
UV reduction method + NED-SAA,
5 ppb
- **ORTOPHOSPHATE (5 min.):**
Molibdate-Antimony, 3 ppb

	NH ₄ -N	NO ₂ -N	NO ₃ -N	PO ₄ -P
Recovery (%)	96	126	94	109
RSD (%)	22.3	1.2	8.6	4.2
LOD (ppb)	6	1	5	3
LOQ (ppb)	20	3	15	10






Analytical procedures for nutrients

	PO ₄ -P [10,11]	NH ₄ -N [9,10]	(NO ₃ + NO ₂)-N [9,10]	NO ₂ -N [9,10]
Filtration	yes	yes	yes	yes
Sampling	90 sec	90 sec	90 sec	90 sec
Sample blank reading	---	yes	---	---
Injection of first reagent	Acid Molybdte	OPA	DTPA	SAA
Sample blank reading	yes	---	yes	yes
UV reduction	---	---	5 min	---
Second reagent injection	Ascorbic acid	---	SAA	NED
Mixing	yes	yes	yes	yes
Third reagent injection	---	---	NED	---
Temperature conditioning	25 °C*	45 °C, 6 min.	---	---
Duration of analysis	5 min	12 min	12 min	5 min
Flow cell	20 mm	10 mm	20 mm	20 mm
Wavelength/ Detection method [9,10,11]	880 nm, Spectrophotometric	370/420nm Fluorometric	525 nm Spectrophotometric	525 nm Spectrophotometric
Range	5-500 µg/L	3-1000 µg/L	10-500 µg/L	2-200 µg/L
Wash cycle	50 sec	50 sec	50 sec	50 sec
Inlet line back wash	---	---	---	5 sec, DIC /acid

* if required under low ambient temperature, water sample temperature could be increased to 25 °C

PO₄

NO₂



The nitrite undergoes a diazotization with fluorosceinamine in an acidic medium. After sufficient incubation the reaction is completed by addition of a strong base. The resulting dye is measured fluorometrically (480nm excitation and 520nm emission) and is directly proportional to concentration¹.

¹ Axelrod, H. *et al*, Anal. Chem., **47**, 922-924 (1975).



MicroLFA: available parameters

- Nutrients: NH_3 , $\text{NO}_3 + \text{NO}_2$, NO_2 and PO_4
- SiO_2
- Total Phosphorus (measurement time 1 h)
- Total Nitrogen (measurement time 30 min.)
- Metal ions like Cr^{6+} , Al, Cu, Iron, Zn, Mn
- Urea.

- **automatic sample blank correction**
 - **automatic washing**
 - **automatic sample dilution allows double scale measurements**
 - **plug-in multi-hydraulic connector available for easy reagents changeover**
 - **compactness and modularity allow easy integration in Ferrybox**
 - **RS-232 protocol compatibility with Micromac-1000 and sondes**
 - **compact dimensions: 270 (H) x 150 (L)x 175 (W) mm, hydraulics / electronics.**
-




Proposed field test under JERICO TransNational Access program



Submitted a proposal to test the MicroLFA modules in the following HZG's facilities:

- **Cuxhaven fixed monitoring station at the Elbe river mouth**
- **Ferrybox system on a regular route.**

0.2 microns filtration unit with auto back-wash





Requirements for a chemical analyzer in a Ferrybox system

- ✓ Long term unattended autonomy
- ✓ Compactness
- ✓ Low reagents and sample consumption
- ✓ Low life-cycle cost
- ✓ Excellent reliability
- ✓ Easy interface with data-logger
- ✓ Low maintenance by non expert users

DONE !



THANK YOU
FOR YOUR KIND ATTENTION